



DISCOVER THE POWER OF LATERAL FLOW TESTS: FAST ALLERGEN DETECTION WITH **INVIRAPID®**

The InviRapid® Lateral Flow Assays allow rapid on-site detection of a variety of food allergens within minutes in a broad range of food matrices, rinse waters and surface swabs.

LATERAL FLOW TESTS HAVE NUMEROUS ADVANTAGES

- Fast results
- Suitable for on-site analysis
- Highly specific
- Validated
- No laboratory equipment needed
- Free from any hazardous substances
- Easy to use, even by inexperienced personnel



The interpretation of results is achieved through visual analysis, utilizing specific antibodies labeled with colloidal gold. For convenient on-site testing, the kits provide all needed reagents and consumables such as test strips, sampling tubes and devices, extraction buffer, and additional accessories.

The InviRapid® test strips can be seamlessly integrated into allergen sanitation and management plans, serving as a valuable component of a Hazard Analysis and Critical Control Points (HACCP) program for all food-related businesses.





GET READY TO **SIMPLIFY YOUR WORKFLOW**
WITH INVITEK DIAGNOSTICS

Detection kit	Catalogue No.
InviRapid® Casein Allergen Ag Test	6022001102
InviRapid® Shellfish Allergen Ag Test	6022002102
InviRapid® Gluten Allergen Ag Test	6022003102
InviRapid® Egg Allergen Ag Test	6022004102
InviRapid® Fish Allergen Ag Test	6022005102
InviRapid® Soya Allergen Ag Test	6022006102
InviRapid® Almond Allergen Ag Test	6022007102
InviRapid® Peanut Allergen Ag Test	6022008102
InviRapid® Sesame Allergen Ag Test	6022009102

CONTACT US

info@invitek.com

INVITEK DIAGNOSTICS PORTUGAL

ALS Life Sciences, SA
Zona Industrial de Tondela, ZIM II, Lote 6
3460-070 Tondela, Portugal

+351 232 817 817


INVITEK DIAGNOSTICS GERMANY

Invitek Molecular GmbH
Robert-Roessle-Str. 10
D-13125 Berlin, Germany

+49 30 9489 2908

Standard package size 20 tests

BE CURIOUS
AND SEE WHAT'S IN FOR YOU



invitek.com