



## DISCOVER THE POWER OF LATERAL FLOW TESTS: FAST ALLERGEN DETECTION WITH **INVIRAPID®**

The InviRapid® Lateral Flow Assays allow rapid on-site detection of a variety of food allergens within minutes in a broad range of food matrices, rinse waters and surface swabs.

### **LATERAL FLOW TESTS HAVE NUMEROUS ADVANTAGES**

- Fast results
- Suitable for on-site analysis
- Highly specific
- Validated
- No laboratory equipment needed
- Free from any hazardous substances
- Easy to use, even by inexperienced personnel



The interpretation of results is achieved through visual analysis, utilizing specific antibodies labeled with colloidal gold. For convenient on-site testing, the kits provide all needed reagents and consumables such as test strips, sampling tubes and devices, extraction buffer, and additional accessories.

The InviRapid® test strips can be seamlessly integrated into allergen sanitation and management plans, serving as a valuable component of a Hazard Analysis and Critical Control Points (HACCP) program for all food-related businesses.



GET READY TO **SIMPLIFY YOUR WORKFLOW** WITH INVITEK DIAGNOSTICS



GET READY TO **SIMPLIFY YOUR WORKFLOW**  
WITH INVITEK DIAGNOSTICS

Detection kit	Catalogue No.
InvitRapid® Casein Allergen Ag Test	6022001102
InvitRapid® Shellfish Allergen Ag Test	6022002102
InvitRapid® Gluten Allergen Ag Test	6022003102
InvitRapid® Egg Allergen Ag Test	6022004102
InvitRapid® Fish Allergen Ag Test	6022005102
InvitRapid® Soya Allergen Ag Test	6022006102
InvitRapid® Almond Allergen Ag Test	6022007102
InvitRapid® Peanut Allergen Ag Test	6022008102
InvitRapid® Sesame Allergen Ag Test	6022009102

Standard package size 20 tests

Contact us

[info@invitek.com](mailto:info@invitek.com)


**INVITEK DIAGNOSTICS PORTUGAL**

Zona Industrial de Tondela, ZIM II, Lote 6  
3460-070 Tondela, Portugal

**INVITEK DIAGNOSTICS GERMANY**

Robert-Roessle-Str. 10  
D-13125 Berlin, Germany

BE CURIOUS  
AND SEE WHAT'S IN FOR YOU



[invitek.com](https://invitek.com)