




INVISCREEN® PCR KITS **NEW USER-FRIENDLY PROTOCOLS**

OUR INVISCREEN® PCR KITS FOR FOOD TESTING AND PATHOGEN DETECTION GIVE YOU FAST AND ACCURATE RESULTS

- New optimized PCR reaction setup - easily pipetted in just 2 steps
- Proven accuracy and sensitivity
- For food allergen testing, GMO detection, foodborne pathogen detection, animal health and infectious diseases testing

EXPERIENCE EFFORTLESS CONVENIENCE

With just 2 easy steps, you can pipette the PCR preparations effortlessly.



Maximize your efficiency in the lab and save valuable time!

Each kit comprises a qPCR master mix and a primer/probe mix. Simply combine both components at 10 µl per reaction, and they are instantly ready to use. No need to add water or spend time preparing complex reaction mixtures.

You can rely on the reliability, accuracy and sensitivity of our InviScreen® kits to give you consistent and reliable results.

BE CURIOUS
AND SEE WHAT'S IN FOR YOU



[invitek.com](https://www.invitek.com)



GET READY TO
SIMPLIFY YOUR WORKFLOW
WITH INVITEK DIAGNOSTICS

SELECTION OF INVISCREEN® KITS

ALLERGEN DETECTION

Detection kit	Catalogue No.
InviScreen® Almond DNA Detection Kit	6012001200
InviScreen® Brazil Nut DNA Detection Kit	6012002200
InviScreen® Cashew DNA Detection Kit	6012003200
InviScreen® Celery DNA Detection Kit	6012004200
InviScreen® Crustacean DNA Detection Kit	6012005200
InviScreen® Hazelnut DNA Detection Kit	6012006200
InviScreen® Macadamia nut DNA Detection Kit	6012007200
InviScreen® Mustard DNA Detection Kit	6012008200
InviScreen® Peanut DNA Detection Kit	6012009200
InviScreen® Pecan nut DNA Detection Kit	6012010200
InviScreen® Pistachio DNA Detection Kit	6012011200
InviScreen® Soya DNA Detection Kit	6012012200
InviScreen® Walnut DNA Detection Kit	6012013200
InviScreen® Fish DNA Detection Kit	6012014200
InviScreen® Coconut DNA Detection Kit	6012015200
InviScreen® Lupin DNA Detection Kit	6012016200

Package size 100 reactions

GMO DETECTION

Detection kit	Catalogue No.
InviScreen® GMO Detection Kit (35S/NOS/FMV)	6014001201
InviScreen® Endogenous Plant Detection Kit	6014002201
InviScreen® CaMV Detection Kit	6014003201
InviScreen® GMO Detection Kit (35S/NOS/FMV/pat/bar)	6014004203

Package size 50 reactions/target

HALAL TESTING

Detection kit	Catalogue No.
InviScreen® Swine Halal Speciation Kit	6013009200

Package size 100 reactions

FOODBORNE PATHOGEN DETECTION

Detection kit	Catalogue No.
InviScreen® Norovirus Detection Kit	6015002200
InviScreen® STEC Detection Kit	6015001200
InviScreen® STEC Serotyping Kit	6015004204 *
InviScreen® O157:H7 Detection Kit	6015003201**
InviScreen® C. botulinum Detection Kit	6015005200

Package size 100 reactions | * 6 x 25 reactions | ** 50 reactions

ANIMAL HEALTH

Detection kit	Catalogue No.
InviScreen® Infectious Brusal Disease Virus Detection kit	6016001200
InviScreen® Infectious Bronchitis Virus Detection kit	6016002200
InviScreen® Mycoplasma Detection Kit	6016003200

Package size 100 reactions

CONTACT US

info@invitek.com

INVITEK DIAGNOSTICS PORTUGAL

ALS Life Sciences, SA
Zona Industrial de Tondela, ZIM II, Lote 6
3460-070 Tondela, Portugal

+351 232 817 817

INVITEK DIAGNOSTICS GERMANY

Invitek Molecular GmbH
Robert-Roessle-Str. 10
D-13125 Berlin, Germany

+49 30 9489 2908